

Premium Restorations

	Premium PFM	Lava™ Classic	Pressed e.max™	e.max CAD™	Crystalite	Suprinity	TriLuxe
Material	Semi-precious or precious alloy coping with layered porcelain overlay	Premium 3M zirconia coping with multi-layered porcelain overlay	Ivoclar™ e.max lithium disilicate	Ivoclar™ e.max lithium disilicate	Monolithic super translucent zirconia	Zirconia reinforced lithium silicate	Feldspar ceramic
Method of manufacture	Traditionally-cast & built up by senior technician	Digitally designed and manufactured then built up with porcelain by hand	Waxed & pressed by hand	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed and manufactured
Aesthetics	★★★★★	★★★★★★	★★★★★★	★★★★★	★★★★★	★★★★★	★★★★★★
Shades / colouring <small>(using Vita Classical, 3D-Master or Chromascop guides)</small>	Internal effects customised shading multi-layered	Internal effects customised shading multi-layered	Customised shading	Customised shading	No internal effects standard shades	No internal effects standard shades	No internal effects standard shades
Strength	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
Add-on options <small>(fit to denture, attachments etc. Charges apply)</small>	✓	Fit to denture, rests and ledges	✗	✗	Fit to denture	✗	Rests and ledges
Benefits	Traditional method choice of alloys 10x microscope fit multi-layer build	Premium aesthetics excellent fit metal-free	Premium aesthetics metal-free	Premium aesthetics excellent fit metal-free	High durability extreme strength metal-free, exceptional accuracy	Perfect mix of aesthetics & strength	High aesthetics, minimal wear
Recommended use	Full mouth, crowns & bridges	Full mouth, crowns & bridges	Full mouth, crowns, inlays/onlays veneers, 3/4 crowns	Full mouth, crowns, inlays/onlays veneers, 3/4 crowns	Full mouth, single crowns	Full mouth, single crowns	Anterior crowns, veneers
Warranty	6 years	1 year	1 year	5 years	10 years	5 years	5 years

Standard Restorations

	Standard PFM	Full Gold Crown	Opalite	Zirconia (Layered)	Lava™ Ultimate	Lava™ Plus	Celtra® Duo	Enamic
Material	Semi-precious or precious alloy coping with porcelain overlay	Non-precious (gold-coloured) or precious (gold) alloy CoCr (silver coloured)	Monolithic translucent zirconia	Zirconia (layered) coping with porcelain overlay	Nano-composite	Monolithic translucent zirconia	Zirconia-Reinforced Lithium Silicate	Hybrid Ceramic
Method of manufacture	Traditionally-cast & built up	Traditionally-cast	Digitally designed and manufactured	Digitally designed & manufactured (Framework) & built up by hand	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed and manufactured
Aesthetics	★★★★☆	★☆☆☆☆	★★★★☆	★★★★☆	★★★★☆	★★★★☆	★★★★☆	★★★★☆
Shades / colouring <small>(using Vita Classical, 3D-Master or Chromascop guides)</small>	No internal effects standard shades	Not required	No internal effects standard shades	No internal effects standard shades	No internal effects standard shades	No internal effects standard shades	No internal effects standard shades	No internal effects standard shades
Strength	★★★★☆	★★★★★	★★★★★	★★★★☆	★★★★☆	★★★★★	★★★★☆	★★★★☆
Add-on options <small>(fit to denture, attachments etc. Charges apply)</small>	✓	✓	Fit to denture	Fit to denture, rests and ledges	✗	Fit to denture	✗	✗
Benefits	Traditional method choice of alloys 5x microscope fit	5x microscope fit exceptional strength	High durability extreme strength metal-free, exceptional accuracy	Good fit metal-free	Accurate adjustability, less abrasion	High durability extreme strength metal-free, exceptional accuracy	High durability extreme strength metal-free, exceptional accuracy	Minimal invasive. shock absorbing
Recommended use	Full mouth, crowns & bridges	Posterior teeth	Full mouth, crowns & bridges	Full mouth, crowns & bridges	Onlays / Inlays	Full mouth, crowns & bridges	Single crowns, inlays, onlays	Full mouth, crowns, inlays, onlays, implants
Warranty	6 years	6 years	15 years	1 year	10 years	15 years	10 years	5 years